

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL SHEET

PRODUCT NAME	ECONOMICAL FACE SHIELD WITH 6 VISORS
CODE	R5
COMPONENTS	PROTECTIVE VISOR: POLYCARBONATE ANTI-FOG VISOR THICKNESS: 0.178
THERMAL CHARACTERISTICS	THE VISOR IS REUSABLE AND NOT STERILISABLE
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKING	PACKING, IN CARTON BOX: 6 PCS THE VISOR COVERS NOSE MOUTH AND EYES THE PROTECTIVE VISOR HAS A DOUBLE PROTECTIVE COATING THE VISOR IS VERY TRANSPARENT
CHEMICAL CHARACTERISTICS	IT RESISTS TO CHEMICAL DISINFECTANTS
SAFETY DATA SHEET	N/A
PRECAUTION	WARNING: USE ONLY AS FACE SHIELD AGAINST SPLASHING LIQUID, WATER OR SALIVA
INSTRUCTIONS FOR USE	SUPPLIED IN "NON-STERILE" CONDITION. AUTOCLAVABLE. THE PRODUCT IS FOR SINGLE USE. AVAILABLE INSTRUCTIONS FOR FITTING INSIDE CARTON BOX.
ANALYSIS RISK REF.D.L.46 24/2/97	AVAILABLE ONLY IN ITALIAN
CLASS CEE EN 166	IS AN INDIVIDUAL PROTECTION DEVICE (DPI)

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	MASCHERA BIANCA ECONOMICA CON 6 VISIERE. ALTAMENTE PROTETTIVA
NS CODICE	R5
MATERIALE	VISIERA: POLICARBONATO ANTI APPANNANTE SPESSORE VISIERA: 0.178 mm
PROPRIETÀ TERMICHE	LA VISIERA E' RIUTILIZZABILE E NON STERILIZZABILE
PROPRIETÀ MECCANICHE, FISICHE, CARATTERISTICHE, CONFEZIONAMENTO	CONFEZIONE: 6 PEZZI LA VISIERA CHE COPRE NASO E BOCCA, HA UN DOPPIO RIVESTIMENTO PROTETTIVO ED E' DOTATA DI UNA ECCELLENTE TRASPARENZA
PROPRIETÀ CHIMICHE	RESISTE A TUTTI GLI AGENTI CHIMICI DISINFETTANTI
SCHEDA DI SICUREZZA	N/A
PRECAUZIONI	ATTENZIONE: LA VISIERA PROTETTIVA PROTEGGE SOLO ED ESCLUSIVAMENTE DA GOCCE E SPRUZZI D'ACQUA O SALIVA
ISTRUZIONI PER L'UTILIZZO	DESTINATA AD ESSERE UTILIZZATA ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO. VIENE FORNITA IN CONDIZIONE "NON STERILE" ED IL PRODOTTO È MONOUSO
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	SI
CLASSIFICAZIONE CEE EN 166	E' UN DISPOSITIVO DI PROTEZIONE INDIVIDUALE (DPI)

N/A = NON APPLICABILE N/N = NON NECESSARIA