

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

TECHNICAL DATA SHEET

PRODUCT NAME	COMPOSTABLE PLASTIC ORTHOBOX
OUR ART. CODE	BIO-BOX
COMPONENTS	BIO-BASED GREEN PE (RECYCLABLE PET)
THERMAL CHARACTERISTICS	
PHYSICAL CHARACTERISTICS	RECYCLABLE 8x7 h 4 cm COLOUR WHITE PACKED IN 10 PCS
INSTRUCTIONS FOR USE	N/A
PRECAUTIONS	DO NOT STORE IN SUNLIGHT OR IN TOO HOT AN ENVIRONMENT

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it

LARIDENTMADE IN ITALY
SINCE 1985

SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	ORTHOBOX MATERIALE PLASTICO COMPOSTABILE
NS CODICE	BIO-BOX
MATERIALE	BIO-BASED GREEN PE (PET RICICLABILE)
PROPRIETÀ TERMICHE	
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	RICICLABILE 8x7 h 4 cm COLORE BIANCO CONFEZIONE 10 PZ
ISTRUZIONI PER L' UTILIZZO	N/A
PRECAUZIONI	NON IMMAGAZZINARE AL SOLE O IN AMBIENTI TROPPO CALDI