

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL DATA SHEET


PRODUCT NAME	WHITE BIODEGRADABLE DISPOSABLE DAPPEN
OUR ART. CODE	BIO-DAPPEN
COMPONENTS	BIO-BASED GREEN PE (RECYCLABLE PET)
THERMAL CHARACTERISTICS	THE PRODUCT CAN'T RESIST AT HIGH TEMPERATURES ABOVE 40°C THE PRODUCT MODIFIES
PHYSICAL CHARACTERISTICS	RECYCLABLE APPLICABLE TO TRAYS COLOUR NATURAL WHITE PACK OF 200 PCS
INSTRUCTIONS FOR USE	ONLY TO BE USED BY A SPECIALIST DOCTOR OR HIGHLY QUALIFIED STAFF DISPOSABLE PRODUCT
PRECAUTIONS	DO NOT STORE IN SUNLIGHT OR IN TOO HOT AN ENVIRONMENT

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	DAPPEN MONOUSO BIODEGRADABILE
NS CODICE	BIO-DAPPEN
MATERIALE	BIO-BASED GREEN PE (PET RICICLABILE)
PROPRIETÀ TERMICHE	PLASTICA BIO-BASED PROVENIENTE DA FONTI RINNOVABILI
PROPRIETÀ MECCANICHE, FISICHE E CARATTERISTICHE	RICICLABILE APPLICABILE AI VASSOI COLORE BIANCO NATURALE CONFEZIONE 200 PZ
ISTRUZIONI PER L' UTILIZZO	DESTINATI AD ESSERE UTILIZZATI ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO IL PRODOTTO E' MONOUSO
PRECAUZIONI	NON IMMAGAZZINARE AL SOLE O IN AMBIENTI TROPPO CALDI